

A flowable seed treatment fungicide for the control of seed borne diseases

A flowable seed treatment fungicide with dual modes of action to give effective control of seed borne diseases of wheat, barley, oats, rye and triticale.

Application

- Anchor can be applied through most seed treatment machines, both batch and continuous types, where secondary mixing is provided
- Machines should be calibrated to apply the recommended rate and this should be distributed evenly between the seeds with good coverage to provide optimum disease control
- Co-application of an equal volume of water with Anchor is likely to improve the evenness of coverage of treated seed
- After application, machinery and pipelines should be washed with clean water
- For the latest information on compatibility (co-application) with other seed treatments contact Certis Europe UK
- Flow of Anchor treated seed through seed drills may be affected relative to untreated seed
- Always recalibrate seed drills before sowing treated seed

Seed quality and storage

- Anchor should be applied to good quality seed with a moisture content less than 16%
- Avoid treating cracked, split or sprouted seed and seed of low germination capacity and vigour
- Seed treated with Anchor should not be retained from one season to the next

Bunt



Formulation	200g/l carboxin and 200g/l thiram		
Crop and max. application rate (Litres per tonne of seed)	Wheat		3.0L
	Winter barley		3.0L
	Spring Barley		3.0L
	Oats		3.0L
	Rye, Triticale		3.0L
Max application rate	300ml/10kg seed	Control	Seed borne diseases
Pack	20 Litre, 200 Litre and 1000 Litre HDPE containers	MAPP no.	16488

USE PLANT PROTECTION PRODUCTS SAFELY. ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE.

For further information with regard to the warning phrases and symbols for this product please refer to the product label. Anchor[™] is a registered trademark of Chemtura Corporation. Anchor[™] contains carboxin and thiram. MAPP no. 16488
Marketed by: Certis 1 Riverside, Suite 5, Granta Park, Great Abington, Cambridge CB21 6AD
Tel : 0845 373 0305 • Fax: 01223 891210



scan here to download our free
Product guide app

or visit m.certiseurope.co.uk from
your smart phone or tablet